



medifa 2000

Multifunctional
Tables



we care.

medifa 2000

- Ergonomic and modern design
- Optimal foot space because of Y-base construction
- Optimal cleaning and disinfection
- User-friendly foot switch conception
- Safety engines against wedge and squeezing hazard
- Electric adjustment of head part as standard function
- Single or parallel movement for height adjustment columns
- Anti- / Trendelenburg position available
- Accessible for handicapped and elderly people
- Various colouring possible



Functionality and Design



Height adjustment from
580 up to 980 mm



Trendelenburg adjustment from
-22° up to +22°



Head section adjustment
from -27° up to +40°



Back section adjustment from
0° up to +46°



Multifunctional table with two-part lying surface

Art.-no. 112802

- Electromechanical height adjustment from 580 to 980 mm
- Electromechanical head section adjustment from -27° up to +40°
- Electromechanical Trendelenburg adjustment from -22° up to +22°
- Frame colour: White (RAL 9010)
- Easy adjustment of powered functions by foot switch
- Upholstery thickness: 70 mm
- Safety load: 225 kg
- Lying surface: 1960 x 650 mm
- Floor space required: 1750 x 620 mm
- Power supply: 230 V / 50-60 Hz

Shown Accessories:

85200

Side rails

83201

Nose slit incl. filling pad



Multifunctional table with three-part lying surface

Art.-no. 112805

- Electromechanical height adjustment from 580 to 980 mm
- Electromechanical head section adjustment from -27° up to +40°
- Electromechanical Trendelenburg adjustment from -22° up to +22°
- Electromechanical back section adjustment from 0° up to +46°
- Frame colour: White (RAL 9010)
- Easy adjustment of powered functions by foot switch
- Upholstery thickness: 70 mm
- Safety load: 225 kg
- Lying surface: 1960 x 650 mm
- Floor space required: 1750 x 620 mm
- Power supply: 230 V / 50-60 Hz

Shown Accessories:

112898

Mobility white

85200

Side rails

81700

Leg holders according to Goepel Classic

81710







Artificial leather cover of leg holders in colour of upholstery

1128530

Hand grips white

28531

Artificial leather covering of hand grips












	Range:	112802	112805
 Height adjustment	580 - 980 mm	●	●
 Trendelenburg	-22° / +22°	●	●
 Head section adjustment	-27° / +40°	●	●
 Back section adjustment	0° / +46°		●
 Mobility		◆	◆
 Safety load		225 kg	225 kg

● Electromechanical ◆ Optional


Technical Data
MPG - Class: I
Protection Class: I
Protection type: IPX4
Power supply: 230 V

- Upholstery**
- Ergonomic implementation of cushion parts
 - Rounded rough edges
 - High quality artificial leather
 - Seamlessly manufactured
 - Optimal possibilities to clean and disinfect
 - Various colourings without any extra charge
 - Upholstery spacers against wedge and squeezing hazard




 -50	 -56	 -63	 -67	 -75
White	Beige	Lime	Lemonade	Agave
 -73	 -61	 -65	 -53	 -58
Apricot	Carnelian	Cherry	Pearl	Stone
 -55	 -62	 -74	 -80	 -81
Dolphin	Light Blue	Atoll	Royal	Purple
 -107	 -83	 -96		
Graphite grey	Black	Black, electrically conductive		

Models 112802



Models 112805



Accessories



Art.-no. 28013

Neck upholstery for attachment at back rest (model 112805)



Art.-no. 83201

Nose slit in head section, incl. filling cushion



Art.-no. 85200

Side rails for attachment of accessories, l= 600 mm, pair



Art.-no. 112898

Mobility, white RAL 9010, with four twin castors



Art.-no. 81001

Support for paper roll holder stick (in combination with 81524)



Art.-no. 81524

Paper roll holder stick, l=650 mm



Art.-no. 81800

Infusion holder, extendable, with 4 hooks and dropping glass, incl. clamp



Art.-no. 1128530

Hand grips, for attachment at pelvis section (model 112805), pair

Art.-no. 28531

Covering of hand grips with artificial leather in upholstery colour, pair

Art.-no. 112899_P

Potential equalisation for connecting protective conductor

Art.-no. 112899_S


Key lock for blocking all electric functions

Accessories



Art.-no. 81700

Leg holders according to Goepel Classic with black foamed PUR knee bolsters, pair, without clamps (61511)



Art.-no. 81710

Artificial leather cover of leg holders in colour of upholstery, pair



Art.-no. 61511

Clamp for fixation of accessories with round rod (18 mm), rotatable



Art.-no. 61211

Arm rest, rotatable and tiltable by means of ball joint, height adjustable, horizontally slidable, incl. clamp and fixture strap, upholstery made of PUR (H = 40 mm), L x W: 400 x 115 mm

medifa

The medifa healthcare group unites competence concerning operating rooms under one roof.

medifa is an owner-run and internationally active complete solution provider in the medical technology.

furniture. In addition, one of the largest German OEM suppliers in the medical technology sector also belongs to medifa.

take place exclusively in Germany, and with over 90% added value in the group of companies, the group guarantees fast reaction times.

The product and service portfolio covers all operating room requirements from modular wall, door, ceiling and cabinet systems (RooSy® Room Systems) for the turnkey operating room, as well as mobile operating tables, operating table accessories and operating lights, right up to examination chairs and medical


medifa GmbH & Co. KG (Finnentrop), medifa hygienic rooms GmbH (Ötigheim) and medifa metall und medizintechnik GmbH in Rastatt are all part of this group of companies. The medifa healthcare group employs approximately 350 employees at three locations within Germany. The planning, engineering and production

The slogan “we care” means for us that we consistently orientate ourselves according to the requirements of our customers. Customers all over the world can rely on top-class consultation, processing and quality from a single source, made in Germany according to certified national and international quality standards.

Products

- Modular operating room systems (RooSy)
- Mobile operating tables and accessories
- Operating lights (H!Light)
- Examination and treatment chairs
- Medical furniture and trolleys
- Customized medical solutions

medifa worldwide



we care.



Version 02-21-EN

Subject to modifications due to technical development.
This document is intended to provide information to an
international audience outside of the US.

www.medifa.com



medifa GmbH & Co. KG

Industriestraße 5
57413 Finnentrop
Germany
fon +49 2721 71 77 -0
fax +49 2721 71 77 -255
info@medifa.com

we care.